Filterless Centifuge Extractor ZI 4016

Standards:

Specification:

- The equipment is in accordance with EN 12697-1, EN 13108
- Centrifuge is designed for proper separation of the filter quickly without fillers or binders containing suspended sediment mixtures
- Is not required no filter, no dispersion of the material, which ensures maximum accuracy so
- Solution is poured into a funnel and falls in container dia test. 70x200mm that rotates, due to centrifugal effect, the liquidrises vertically, leaving the pot filler and mineral particles
- The centrifuge is supplied complete with aluminum tank, two sites of 2mm and 0.063mm
- Rotation speed is 11,500 rpm and extractive capacity is up to 100g of filler per test
- 230v, 50hz, 600w

